


Helpful Hints

Congratulations on your purchase of the HTP-1500. Before we begin, there are a few things you should know about this system and heat therapy in general. ►




Fill



Connect



Apply



On



Set Temp

- 1 We recommend that you plug the HTP-1500 into a surge protector if possible. It's not required, but it could help protect the unit during an electrical storm.
- 2 After you unpack the unit and remove the cap, you may see some water inside the jug. That's okay, the unit was tested at the factory and not all the water could be removed. Your HTP-1500 is brand new!
- 3 Before you attach the pad to the hose, be sure to uncoil the hose from the HTP-1500. Operating instructions are located on the pump's housing.
- 4 Do not run the HTP-1500 without unwrapping the hose from its storage position.
- 5 You should hear a "click" when you attach the pad to the hose. If you do not hear a "click" disconnect it and try again. That "click" indicates you have a good connection and water will flow properly.
- 6 The HTP-1500 is preset to 105°F. This is a safe temperature. However, your doctor may advise that you lower the temperature if you will be using the pad for long periods of time (greater than 30 minutes).
- 7 Moisturizers are okay, but **NEVER COMBINE** Ben Gay®, Icy Hot®, Flexall 454® or other heating creams (analgesics) with your use of the HTP-1500; injury to your skin could result.
- 8 When using Heat Therapy, your doctor may recommend checking your skin condition at least every 20 minutes.
- 9 If you notice redness or experience discomfort, remove the pad and turn the HTP-1500 off. If you want to use the HTP-1500 again, consult your doctor, he may recommend that you lower the temperature or discontinue the heat therapy.

Adroit
Medical Systems
adroitmedical.com



Be sure to use your HTP-1500 properly and with caution. Never combine heat therapies. Consult your doctor and discontinue use if you experience discomfort or unusual redness in your skin.