

Hypothermia Blankets Heat Therapy Pads

Thermoregulation & Temperature Therapy Products

adroitmedical.com

Hypothermia Blankets

Hypo-hyperthermia blankets have never been better! Made of polyurethane, these single-patient use blankets are soft to the touch and conform very well to irregular shapes. They are also lightweight and will remain soft even at low temperatures. Fold them or wrap them and don't worry about water not getting through, the unique pinwheel design allows water to flow in hundreds of directions; ensuring good performance.



Super Soft
Ultra Thin
Lightweight
Fold-able
Wrap-able

Soft-Temp® blankets and pads are intended for use by or on the order of a physician only. Applications may include; prevention of inadvertent hypothermia in the surgical patient, re-warming patients after surgery, fever reduction, induced hypothermia and therapeutic hypothermia.



Madison Gammons, grandson of original hypo-hyperthermia blanket inventor, Clifford Gammons.
U.S. Patent No. 4,149,541

Heat Therapy Pads

Soft-Temp® heat therapy pads are so nice that patients have actually called from their hospital bed to ask "Where can I buy one of these?" And for good reason - Effective heat therapy can be a good source of pain relief without the side effects of medications. Soft-Temp® pads can be folded, wrapped, and or wetted on the white side to provide moist heat therapy.



Small / Infant



Medium



Large

Hypothermia Blankets Heat Therapy Pads

Product Information

Soft-Temp® blankets and pads are made from polyurethane; a soft, rubber-like plastic that will not make or leave indentations on the skin. Our clever pin-wheel design ensures unrestricted water flow and an even dispersion of the patient weight; thus reducing the likelihood of irritations at bony prominences. The internal pressures of our blankets and pads are also very low (usually 3 to 5 psi). The unique design and soft materials allow our blankets to contour easily to irregular surfaces and literally drape over a patient's body (you won't see that with polyethylene or vinyl-based products). Soft-Temp® blankets are also latex free and safe to use in direct contact with the patient's skin. Gel mats may be used, but we don't recommend putting them between the patient and the blanket as it could hamper heat-transfer.

Blankets & Pads

- ST-060** Adult blanket, fabric covering, 26" x 66"
- ST-030** Pediatric blanket, fabric covering, 26" x 33"
- ST-024** Large pad, fabric covering, 20" x 26"
- ST-020** Medium pad, fabric covering, 16" x 22"
- ST-017** Small pad, fabric covering, 13.5" x 18.5"
- ST-003** Neck collar, fabric covering, 4" x 19"*

*Other blankets and pads are available, see web site for more information. All dimensions are approximate and before filling. Blankets and pads shrink slightly when filled. *ST-003 not shown.*

Heat therapy is used to treat common pain disorders, including low back pain, pain from arthritis, trauma, muscle strain, wrist pain, and menstrual cramping. Continuous low level heat therapy (approx. 104°F / 40°C) is clinically proven to alleviate pain, relax muscles and increase tissue blood flow.¹ Both independent research and the FDA/CPSC recommend water-circulating systems over electric heating pads for consistent heating and better patient safety.²

1 American Pain Society, Nov 2004 2 FDA Public Health Advisory, 1996

Device Advise

As with any medical device, hypo-hyperthermia blankets and heat therapy pads should be used within the parameter of clinical judgment. Avoid fluids gathering beneath the patient and check skin condition frequently. Patients at greatest risk are those unconscious, with poor circulation, or pre-existing skin conditions.

