

Foot Wrap

Item# ST-019



Adroit
Medical Systems
800-267-6077

Thermoregulation & Temperature Therapy Products

adroitmedical.com

Soft-Temp® temperature therapy wraps can be used with hot or cold therapy. They deliver 40% more energy than ordinary pads and are perfect for post-op swelling, physical therapy and chronic pain applications. The ST-019 Foot Wrap can be used on the left or right foot or elbows.

Infrared photos before and after demonstrates the foot wrap's ability to increase foot circulation. A 154% increase in relative pulse strength was measured after a 60 minute heat therapy treatment. A 45% increase was measured 30 minutes after heat therapy was removed.



All Soft-Temp Temperature Therapy Pads and Wraps include Velcro fasteners and foam insulation to prevent condensation and maintain temperature. Most are available non-sterile and sterile*. Soft-Temp pads and wraps are intended for use with **HTP-1500** continuous low level heat therapy systems and **Mobile ICE** cold therapy systems. Mobile ICE cold therapy accessories include an automatic on-off switch for time controlled alternating therapy.**

*Call for availability. **Timer available with Mobile ICE model MI-12T only.

ST-019 Product Details: 17"x19", Velcro Fasteners, Hot or Cold, Ambidextrous, Insulated, Quick Disconnects, Hand Washable.

Adroit Medical Systems 800-267-6077 865-458-8600 adroitmedical.com or 4mobileice.com