

Convective Air Warming System

[adroitmedical.com](http://adroitmedical.com)

The Soft-Air OR-110 warming system is easy to use, heats up fast and runs quietly. Trouble free 2-year warranty with extended warranty options available. Widest variety of blankets and accessories available, see reverse side for details.




**HEPA Filtration**  
**UL 60601-2-3rd Ed.**

## Convective Air Warming System


[adroitmedical.com](http://adroitmedical.com)

The Soft-Air warming system uses double-filtration to maximize patient safety. Soft-Air uses a premium HEPA filter at the “blower” unit and a proprietary nonwoven blanket material that filter air a second time at the patient. Soft-Airs double-filtration prevents air jets and turbulence associated with “holes” in a blanket.

Soft-Air has the most comprehensive variety blankets (47) to meet the need of any surgical procedure. A few of the more common blankets are shown below. Soft-Air blankets are available in standard (white-on blue) or reflective (silver on blue).



The Soft-Air warming unit (OR-110 blower) has a complete set of optional accessories, including replacement hose, HEPA filter, Hose-Grip, Hose Hook and Cart with wire basket.



Soft-Air warming blankets are covered or may be covered by one or more patents awarded to Adroit; those include US Patents 6,689, 155; 6,699,270; 6,709,447; 6,827,729; 6,994,720; 7,022, 130; 7,066,949; 7,291, 163; and or additional pending patents.